

**Integra<sup>®</sup> MicroFrance<sup>®</sup>**  
Integra Laparoscopy catalogue

<b>Table of contents</b>	ii-iii	<b>Rotatable handles, for <math>\varnothing</math> 5 and 10 mm non-insulated tube</b>	44
<b>Warranty and Customer Service Information</b>	iv	Angled handle	44
<b>Indications - Contraindications</b>	v	Axial handle	44
<b>MicroFrance® modular Instruments</b>	vi-vii	Ergonomic handle	44
<b>Fenestrated forceps <math>\varnothing</math> 5 mm</b>	2-6	<b>Non-insulated tubes</b>	45
<b>Graduated fenestrated forceps <math>\varnothing</math> 5 mm</b>	7	for rotatable instrument $\varnothing$ 5mm	45
<b>Fenestrated forceps <math>\varnothing</math> 10 mm</b>	8-9	for rotatable instrument $\varnothing$ 10mm	45
<b>Babcock forceps</b>	10-11	<b>Non rotatable handles, <math>\varnothing</math> 5 mm and insulated tubes</b>	46
<b>Alligator forceps <math>\varnothing</math> 5 mm</b>	12-13	Angled handle	46
<b>Teeth forceps</b>	14-16	Axial handle	46
<b>Dissector <math>\varnothing</math> 5 mm</b>	17-19	Black insulated tube for non-rotatable instrument	46
<b>Dissector <math>\varnothing</math> 10 mm</b>	20-21	<b>Non rotatable ergonomic handles and non insulated tubes</b>	47
<b>Needle-Holder 5 mm</b>	22-23	<b>Bipolar cables</b>	48
<b>Needle-Holder 10 mm</b>	24	<b>Non rotatable bipolar forceps, <math>\varnothing</math> 5 mm</b>	49-52
<b>Knot Pusher, Knot guide scissors</b>	25	<b>Monopolar hook</b>	53
<b>Scissors <math>\varnothing</math> 5 mm</b>	26-28	<b>Retractable knives, Palpation probe</b>	54
Black sheath for metal tube	27	Knife blade	54
Complete non modular instruments	28	Monopolar cable	54
Insulated scissors insert	27	Monopolar palpation probe	54
Metal tube	27	Retractable knife	54
Modular instrument	26	Retractable monopolar electrode	54
Wire insert	27	<b>Retractors - Palpation probes</b>	55
<b>Scissors <math>\varnothing</math> 10 mm</b>	29	Articulated retractor	55
<b>Vascular clamp <math>\varnothing</math> 5 mm</b>	30	Flat palpation probe	55
<b>Vascular clamp <math>\varnothing</math> 12 mm</b>	31	Flexible retractor	55
<b>Vascular clamp <math>\varnothing</math> 10 mm</b>	32-33	Round palpation probe	55
<b>Intestinal clamp</b>	34-35	<b>Flexible tension arm</b>	56
<b>Aortic percutaneous clamp</b>	36-37	Instrument holder	56
<b>Temporary Clamp, Clamp applicator</b>	38	Operating table clamp	56
<b>Clips applier, Circular suturing holder insert</b>	39	<b>Electric irrigation suction pump</b>	57
<b>Rotatable angled handles, <math>\varnothing</math> 5 or 3.5 mm</b>	40	Cleaning sleeve for telescope	57
<b>Rotatable axial handles, <math>\varnothing</math> 5 or 3.5mm</b>	41	Electric handle for pump	57
<b>Rotatable ergonomic handles, <math>\varnothing</math> 5 or 3.5 mm</b>	42	Suction irrigation tubes for handle	57
<b>Tubes, <math>\varnothing</math> 5 mm for rotatable handles, <math>\varnothing</math> 5 or 3.5 mm</b>	43	<b>Suction irrigation tube</b>	58
Insulated tube	43	Handle for modular suction irrigation tube	58
Metal tube and pack of 3 black sheaths	43	Modular suction irrigation tube	58
		Non modular suction irrigation tube	58
		<b>Needles - Wire-pass</b>	59
		<b>Cholangiography - Cholescopy</b>	60

<b>Modular trocars - 3.5 mm</b>	61-62	<b>Mini invasive instrumentation, 5mm diameter, stylet handle</b>	82
Trocars with metal cannula with magnetic valve	62	Curved scissors, titanium jaws	82
Trocars with metal cannula with silicone valve	62	Forceps with teeth	82
<b>Modular trocars - 5.5 mm</b>	63	Needle-holder, titanium jaws	82
<b>Modular trocars - 10.5 mm</b>	64	Needle-holder, tungsten jaws	82
Modular trocar 10,5mm	64	<b>Henry set - 3.0mm diameter</b>	83
Trocar tip palpator probe	64	Fenestrated forceps	83
<b>Modular trocars - 12/13 mm</b>	65	Forceps, alligator jaws	83
<b>Trumpet valve trocars</b>	66	Monopolar hook	83
Spring for Trumpet valve trocars	66	Round probe	83
Trumpet valve trocars	66	Scissors	83
<b>Automatic valve trocars</b>	67	Trocar	83
<b>Trocars accessories</b>	68	Trocar holder ring	83
Double retractor	68	<b>Terris set</b>	84
Reducer sleeve for trocar	68	Double retractor	84
Sealing cap for trocar	68	Elevator	84
<b>Thoracic Instrumentation</b>	69	Nerve suction with finger control	84
Fenestrated forceps	69	Peanut holder	84
Forceps, curette jaws	69	Suction tube elevator	84
Insert, curette jaws	69	<b>Scheye rectal forceps</b>	85
Trocar	69	Axial handle	85
<b>Irrigation bipolar forceps</b>	70	Insert, biopsy jaws and tube 5mm diam	85
<b>3.5mm diameter instrumentation</b>	71-75	<b>Sterilization case, other accessories</b>	86
Adapter cleaning luer-lock	72	Cap	86
Bipolar forceps	74	Stainless steel sterilization case	86
Complete non modular instrument	73	Surgical instrument oil	86
Forceps without port	75	Tray	86
Knot pusher	75	<b>Numerical index</b>	87-91
Modular instrument	72	<b>Table for manufacturer part number equivalent</b>	92
Monopolar hook	75		
Palpation probe	75		
Suction irrigation tube	75		
<b>Pediatric instrumentation, 3.5mm diameter</b>	76-81		
Adapter cleaning luer-lock	80		
Bipolar forceps	79		
Complete non modular instrument	81		
Composite trocar	77		
De Lagausie monopolar instrument	81		
Metallic trocar	76		
Modular instrument	80		
Monopolar hook	78		
Needle-holder	81		
Retractable knife	78		
Suction irrigation tube	78		

## Integra Limited Warranty

Please refer to the warranty included in the packaging of each product to know the applicable warranty.

## Contacts Customer service:

### International

+33 (0) 4 37 47 59 50  
+33 (0) 4 37 47 59 25 (fax)  
csemea@integralife.com

### France

+33 (0) 4 37 47 59 10  
+33 (0) 4 37 47 59 29 (Fax)  
custservfrance@integralife.com

### Belgium & Luxembourg

+32 (0) 2 257 4130  
+32 (0) 2 253 2466 (Fax)  
custsvcbenlux@integralife.com

### United Kingdom

+44 (0) 1264 312 725  
+44 (0) 1264 312 821 (Fax)  
custsvcs.uk@integralife.com

### Netherlands

+31 (0) 852 083 167  
+31 (0) 207 093 627 (Fax)  
custsvcnetherlands@integralife.com

### Italy

+39 (0) 2 577 8921  
+39 (0) 2 575 113 71 (Fax)  
custsvcitaly@integralife.com

## MicroFrance® Laparoscopy Instrumentation

### INDICATIONS

- All the medical devices showed in this catalogue have indications for diagnosis and/or laparoscopy surgery.
- Sizes are not contractual and are given for information purposes only.

### CONTRAINDICATIONS


- Do not use in indications other than the indications mentioned in the DFU or Labeling.

### How to read this catalogue

Title and/or product family      Symbol  
"Complete instrument"

"Modular" symbol

Instrument to a 1/1 scale except where mentioned otherwise

Fenestrated forceps # 5 mm				
Modular			Handle	
Jaws	Working length			
Scale (x)	Short	Standard	Large	
Johann 200mm	CEV8625CU 345 mm	CEV8625CU 345 mm	CEV8625CU 345 mm	 Rotatable 40 -> 45
Concave 200mm	-	CEV8625CU 345 mm	-	
Curved jaws 200mm	Bertrand CEV8625CU 345 mm	Bertrand CEV8625CU 345 mm	Bertrand CEV8625CU 345 mm	 Non rotatable 46 - 47
Johann 200mm	-	CEV8625CU 345 mm "Sintermetal" design	-	
Double-action jaws 200mm	CEV8625DU 345 mm	CEV8625DU 345 mm	CEV8625DU 345 mm	

Fenestrated forceps # 5 mm			
NON Modular			
Jaws	Handle	Working length	
Scale (x)		Short	Standard
Johann 200mm	Monopolar with ratchet	Johann CEV8650 345 mm	-
Johann 200mm	With spring	Coyot CEV8638C 345 mm	-
Johann 200mm	A coimelline	Montiput CEV865 240 mm	-
Johann 200mm	A coimelline	CEV865C 240 mm	CEV865 345 mm
Johann 200mm	With long ratchet	-	CEV866 345 mm
Concave 200mm	Monopolar with ratchet	-	CEV869U 345 mm
Curved jaws 200mm	A coimelline	-	Bertrand CEV863C 345 mm

DOCUMENT FOR USE IN EUROPE, MIDDLE-EAST AND AFRICA ONLY.

Except contrary indications, pages of this catalogue are presented as below :

- left page : modular instruments
- right page : complete non modular instrument

In the majority of cases, modular instruments are linked to complete instruments. If we don't have any correspondence, you will have an «X» on the line.


Example : CEV8625CU doesn't have complete version.

In the modular inserts table, handles to be used are precised with specific pages. Tubing compatible with handles are mentioned on pages close to handles' pages.

Please refer to modular instruments composition on following page.


MicroFrance® laparoscopic dismantled instruments are composed of three pieces : Jaws, tube and handle..

The black portion of the left tab indicates the section featured in the current page :



In the majority of cases, modular instruments are linked to complete instruments


## Handles



**\*\* Grey rotatable ring for metal tube or non insulated tube  $\varnothing$  5mm and 10mm**






























**\* Black rotatable ring for  $\varnothing$  5mm insulated tube**

Jaws 










Modular scissors insert : insulated scissors insert blades and wire insert for scissors insert

Tubes 

Non rotatable instrument	Rotatable instrument	Insulated tube	Non insulated tube	Tube with independant sheath	Used for bipolar instrument
					
					
					
					
					
























<b>Modular</b>				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 <b>Johann</b> 20x5mm	CEV8625T1U 245 mm	CEV625T1U 345 mm	CEV9625T1U 445 mm	 <b>Rotatable</b> 40 -> 45
 Concave 20x5mm	-	CEV670T1U 345 mm	-	<b>Non rotatable</b> 46 - 47
 Curved jaws 20x5mm	<b>Bertrand</b> CEV8625CU 245 mm	<b>Bertrand</b> CEV625CU 345 mm	<b>Bertrand</b> CEV9625CU 445 mm	
 <b>Johann</b> 20x5mm	-	CEV625P1U 345 mm "Switzerland" design	-	
 Double-action jaws 20x5mm	CEV8625DU 245 mm	CEV625DU 345 mm	CEV9625DU 445 mm	



**NON Modular**










**Complete instruments**










Jaws Scale 1/1	Handle	Working length	
		Short	Standard
 <b>Johann</b> 20x5mm	 With ratchet	<b>Johann</b> CEV405U 365 mm	
	 With spring	– <b>Gayet</b> CEV405BG 365 mm	
	 With ratchet	<b>Montupet</b> CEV815 240 mm –	
	 With ratchet	CEV915G 240 mm	CEV525 345 mm
	 With long ratchet	–	CEV725 345 mm
 Concave 20x5mm	 With ratchet	–	CEV470U 365 mm
 Curved jaws 20x5mm	 With ratchet	–	<b>Bertrand</b> CEV525C 345 mm

	<b>Modular</b>				Page #	
	Inserts			Working length		Page #
	Jaws Scale 1/1	Short	Standard			
 15x5mm	-	CEV4o6IU 340 mm	-	    <b>Rotatable</b> <b>40 -&gt; 45</b>  <b>Non rotatable</b> <b>46 - 47</b>    		
 <b>Desmaizieres</b> 20x3mm	CEV84o7IU 245 mm	CEV4o7IU 345 mm	CEV94o7IU 445mm			
 40x5mm	CEV866oT1U 270 mm	CEV66oT1U 370 mm	CEV966oT1U 470 mm			
 Concave 40x5mm	-	CEV671T1U 370 mm	-			
 <b>Crovella</b> 40x5mm smooth jaws	-	CEV661IU 370 mm	-			
 <b>Crovella</b> 30x5mm smooth jaws	-	CEV662IU 335mm	-			







**NON Modular**

**Complete instruments**

Jaws Scale 1/1	Handle	Working length	
		Standard	
 15x5mm	 With ratchet	CEV406U	360mm
 <b>Desmaizieres</b> 20x3mm	 With ratchet	<b>Desmaizieres</b> CEV407U	365mm
 40x5mm	 With ratchet	<b>Desplantez</b> CEV460U	385mm
	 Axial, with long ratchet	CEV760	365mm
 Concave 40x5mm	 With ratchet	CEV470T1U	385mm

<b>Modular</b>				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Round double-action jaws 15x5mm	-	CEV433DU 350 mm	-	   <b>Rotatable</b> 40 -> 45  <b>Non rotatable</b> 46 - 47   
 Smooth jaws 20x5mm	CEV8625EU 245 mm	CEV625EU 345 mm	-	
 <b>Porpiglia</b> one fenestrated jaws 20x5mm	-	CEV625PU 345 mm	-	
 Concave double-action jaws 15x5mm	CEV8670DU 240 mm	CEV670DU 340 mm	CEV9670DU 440 mm	

















<b>Modular</b>		
Inserts		Handle
Jaws Scale 1/1	Working length	Page #
 Round jaws, Double-action jaws	<b>CEV433-1</b> 370 mm	  <b>Rotatable</b> <b>44 -&gt; 45</b>  
 35x10mm	<b>CEV610-1</b> 375 mm	
 <b>Hüscher</b> Round jaws, Double-action jaws	<b>CEV433H</b> 380 mm	
 Round jaws, Double-action jaws with tube CEV864910B (see page 45)	<b>CEV8434I</b> 295 mm Short	

**NON Modular**

**Complete instruments**



Jaws Scale 1/1	Handle	Working length
 Round jaws, Double-action jaws	 Axial, with ratchet	<b>CEV533-2</b> 370 mm
 35x10mm	 Angled, with long ratchet	<b>CEV410</b> 360 mm
	 Axial, with ratchet	<b>CEV530</b> 370 mm
 Duval	 Angled, with long ratchet	<b>CEV432</b> 365 mm
	 Axial, with ratchet	<b>CEV507</b> 375 mm
 Round jaws	 Angled, with long ratchet	<b>CEV433</b> 365 mm








<b>Modular</b>						
Inserts		Handle				
Jaws Scale 1/1	Working length					
		Page #				
 Ø10mm	<b>Delaitre</b> CEV602-1 370 mm					
		 <b>Rotatable</b> <b>44 - 45</b>				
		 <b>Rotatable</b> <b>40 -&gt; 45</b> <b>Non rotatable</b> <b>46 - 47</b>				
		 <b>Delaitre</b> Ø5mm				
		<table border="0"> <tr> <td style="text-align: center;">Short</td> <td style="text-align: center;">Standard</td> </tr> <tr> <td style="text-align: center;">           CEV8621T1U            250 mm         </td> <td style="text-align: center;">           CEV621T1U            350 mm         </td> </tr> </table>	Short	Standard	CEV8621T1U 250 mm	CEV621T1U 350 mm
Short	Standard					
CEV8621T1U 250 mm	CEV621T1U 350 mm					













## NON Modular

### Complete instruments






















Jaws Scale 1/1	Handle	Working length				
 Ø10mm	 Angled, with long ratchet	<b>Delaitre</b> CEV420 360 mm				
	 Axial with ratchet	<b>Delaitre</b> CEV502 370 mm				
 Ø5mm	 Angled, with long ratchet	<table border="0"> <tr> <td>Short</td> <td>Standard</td> </tr> <tr> <td> <b>Montupet</b>                      CEV821G                      245 mm                 </td> <td> <b>Delaitre</b>                      CEV421                      370 mm                 </td> </tr> </table>	Short	Standard	<b>Montupet</b> CEV821G 245 mm	<b>Delaitre</b> CEV421 370 mm
	Short	Standard				
	<b>Montupet</b> CEV821G 245 mm	<b>Delaitre</b> CEV421 370 mm				
 Axial, with ratchet	– <b>Delaitre</b> CEV521 350 mm					
 Axial with long ratchet	– <b>CEV721</b> 350 mm					

<b>Modular</b>				
Inserts			Handle	
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 <b>Manhès</b> 10mm	-	CEV614T1U 335 mm	CEV9614T1U 435 mm	    <b>Rotatable</b> 40 -> 45  <b>Non rotatable</b> 46 - 47   
 <b>Manhès</b> 15mm	CEV8615T1U 240 mm	CEV615T1U 340 mm	CEV9615T1U 440 mm	
 <b>Manhès</b> 20mm	-	CEV620T1U 345 mm	CEV9620T1U 445 mm	
 <b>Dauphin</b> 20x3mm	-	CEV617IU 340 mm	CEV9617IU 440 mm	
 Serrated jaws inserts, 15x5mm	CEV8624IU 240 mm	CEV624IU 340 mm	CEV9624IU 440 mm	

## NON Modular













### Complete instruments





Jaws Scale 1/1	Handle	Working length	
		Standard	
 10mm	 With ratchet	<b>Manhès</b> CM110U 355 mm	
 15mm	 With ratchet	<b>Manhès</b> CM115U 360 mm	
	 Axial, with ratchet	<b>CEV515</b> 340 mm	
 20mm	 With spring	<b>Manhès</b> CEV220U 365 mm	
	 With ratchet	<b>Manhès</b> CM120U 365 mm	

				<b>Modular</b>					
		Inserts			Handle				
Jaws Scale 1/1			Working length			Page #			
 <p style="text-align: center;">Ø10mm Grasping insert</p>					<p><b>Claw</b> CEV606-1 370 mm</p>			 <p><b>Rotatable</b> 44 -&gt; 45</p> 	
		Short	Standard	Large					
 <p style="text-align: center;"><b>Allis</b> Grasping insert Ø5mm</p>		CEV8611T1U 235 mm	<b>Allis</b> CEV611T1U 335 mm	<b>Allis</b> CEV9611T1U 435 mm		   <p><b>Rotatable</b> 40 -&gt; 45</p>			
 <p style="text-align: center;"><b>Manhès</b> Cobra Ø5mm</p>		-	<b>Manhès</b> CEV612T1U 340 mm	-		<p><b>Non rotatable</b> 46 - 47</p> 			
 <p style="text-align: center;"><b>Manhès</b> Tiger Ø5mm</p>		CEV8613T1U 240 mm	<b>Manhès</b> CEV613T1U 340 mm	-					

## NON Modular







### Complete instruments








Jaws Scale 1/1	Handle	Working length	
		Standard	
 Ø10mm	 Angled with spring	<b>Claw</b> CEV321 360 mm	
	 Axial, with ratchet	<b>Claw</b> CEV506 370 mm	
 <b>Allis</b> Grasping forceps Ø5mm	 With spring	<b>Allis</b> CEV211U 355 mm	
	 With ratchet	<b>Allis</b> CM111U 355 mm	
	 Axial, with ratchet	<b>Allis</b> CEV511 335 mm	
 <b>Manhès</b> Cobra Ø5mm	 With ratchet	<b>Manhès</b> CM112U 355 mm	
	 Axial, with ratchet	<b>Manhès</b> CEV512 335 mm	
 <b>Manhès</b> Tiger Ø5mm	 With ratchet	<b>Manhès</b> CM113U 360 mm	

<b>Modular</b>				
Inserts			Handle	
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 <b>cobra</b> Double-action jaws Ø5mm	<b>CEV8612DU</b> 240 mm	<b>CEV612DU</b> 340 mm	<b>CEV9612DU</b> 440 mm	 <b>Rotatable</b> <b>40 -&gt; 45</b>
 <b>Vallancien</b> Fenestrated jaws insert with teeth Ø5mm 20x5mm	<b>CEV8616IU</b> 245 mm	<b>CEV616IU</b> 345 mm	-	<b>Non rotatable</b> <b>46 - 47</b> 



## Modular

Inserts		Handle
Jaws Scale 1/1	Working length	Page #
 <p><b>Delaitre</b> Double-action jaws, 22mm Ø5mm</p>	<p><b>CEV641DU</b> 355 mm</p>	 <p><b>Rotatable</b> 40 -&gt; 45</p> <p><b>Non rotatable</b> 46 - 47</p> 
 <p>Double-action jaws, 20mm Ø5mm</p>	<p><b>CEV641FU</b> 350 mm</p>	
 <p><b>Allal</b> Double-action jaws 15mm Ø5mm</p>	<p><b>CEV642TFU</b> 345 mm</p>	
 <p>Double-action jaws, angled 20mm jaws 90deg Ø 5mm</p>	<p><b>CEV643IU</b> 345 mm</p>	







<b>Modular</b>				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Short curved Double-action jaws, 12mm	-	<b>CEV532IU</b> 335 mm	-	    <b>Rotatable</b> <b>40 -&gt; 45</b>  <b>Non rotatable</b> <b>46 - 47</b>
 Medium curved Double-action jaws, 15mm	<b>CEV8641T1U</b> 245 mm	<b>CEV641T1U</b> 345 mm	<b>CEV9641T1U</b> 445 mm	
 Long curved Double-action jaws, 18mm	<b>CEV8642T1U</b> 245 mm	<b>CEV642T1U</b> 345 mm	<b>CEV9642T1U</b> 445 mm	  









**NON Modular**

**Complete instruments**
















Jaws Scale 1/1	Handle	Working length
 <p>Short curved Double-action jaws, 12mm</p>	 <p>Axial, with ratchet</p>	<p><b>CEV532</b> 340 mm</p>
 <p>Medium curved Double-action jaws, 15mm</p>	 <p>Axial, with ratchet</p>	<p><b>CEV531-2</b> 340 mm</p>
 <p>Long curved Double-action jaws, 18mm</p>	 <p>Axial, with ratchet</p>	<p><b>CEV531T3U</b> 350 mm</p>

<b>Modular</b>		
Inserts		Handle
Jaws Scale 1/1	Working length	Page #
 <p><b>Delaitre</b> Medium curved 50mm</p>	<p><b>CEV601-1</b> 380 mm</p>	
 <p><b>Delaitre</b> Medium curved - Double-action jaws, 50mm</p>	<p><b>CEV651-1</b> 380 mm</p>	  
 <p><b>Pignata</b> Short curved 90deg, 45mm</p>	<p><b>CEV601P</b> 380 mm</p>	<p><b>Rotatable</b> 44 -&gt; 45</p>
 <p><b>Puglisi</b> Long curved Double-action jaws, 60mm</p>	<p><b>CEV652l</b> 380 mm</p>	

## NON Modular

















### Complete instruments








Jaws Scale 1/1	Handle	Working length
 <p>Medium curved 50mm</p>	 <p>Angled</p>	<p><b>Delaitre</b> CEV430 370 mm</p>
 <p><b>Delaitre</b> Medium curved - Double-action jaws, 50mm</p>	 <p>Axial, with ratchet</p>	<p><b>Nouaille</b> CEV501 380 mm</p> <p><b>CEV501-2</b> 380 mm</p>

<b>Modular</b>				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Tungsten straight jaws	-	CEV500TU 335 mm	CEV9500TU 435 mm	   <b>Rotatable</b> <b>40 -&gt; 45</b>  <b>Non rotatable</b> <b>46 - 47</b>  
 Tungsten curved left jaws	CEV8538TU 240 mm	CEV538TU 340 mm	-	
 Tungsten curved right jaws	CEV8544IU 240 mm	CEV544IU 340 mm	CEV9544IU 440 mm	
 Titanium automatic jaws	CEV8650TU 235 mm	CEV650TU 335 mm	CEV9650TU 435 mm	

## NON Modular







### Complete instruments

Jaws Scale 1/1	Handle	Working length		
		Short	Standard	Large
 Tungsten straight jaws	 Axial, with ratchet	CEV965 230 mm	CEV500T5 335 mm	–
	 Large ergonomic handle with ratchet up	–	CEV700T5 335 mm	–
	 Large ergonomic handle with ratchet up	–	CEV700T5D 335 mm	CEV9700T5D 450 mm
 Tungsten curved left jaws	 Axial, with ratchet	–	CEV538T5 335 mm	–
	 Large ergonomic handle with ratchet up	–	CEV738T5 335 mm	CEV9738T5 455 mm
	 Fine ergonomic handle	–	CEV738-5 350 mm	–
	 Fine ergonomic handle	–	CEV738D5 350 mm	–
 Tungsten curved right jaws	 Axial, with ratchet	–	CEV544 335 mm	–
	 Large ergonomic handle with ratchet up	–	CEV744T5 335 mm	CEV9744T5 455 mm
 Titanium automatic jaws	 Axial, with ratchet	–	<b>Mouiel</b> CEV540-5 335 mm	–
	 Large ergonomic handle with ratchet up	–	CEV650T5D 335 mm	–
	 Fine ergonomic handle	–	CEV740-5 350 mm	–








<b>NON Modular</b>			
Complete instruments			
Jaws Scale 1/1	Handle		Working length
 <b>Nouaille</b> Straight jaws	 Axial, with ratchet		<b>Nouaille</b> <b>CEV500-A10</b> 345 mm
 Left curved jaws	 Axial, with ratchet		<b>CEV538-10</b> 345 mm
 Straight automatic jaws	 Axial, with ratchet		<b>Mouiel</b> <b>CEV540-10</b> 345 mm
	 Fine ergonomic handle		<b>CEV640-10</b> 340 mm














**Knot Pusher - 5 mm**















		<b>NON Modular</b>	
		Complete instruments	
	<p><b>Delaitre</b> CEV442 L 340 mm</p>		<p>Working length «L» Scale 1/3</p>
	<p>CEV442E L 300 mm</p>		
	<p>CEV442Z L 350 mm</p>		

**Knot guide scissors - 5 mm**









		<b>NON Modular</b>		Working length			
		Complete instruments					
Jaws Scale 1/1		Handle		Short	Standard	Large	
				-	<p>CEV445-5 355 mm</p>	-	
<p><b>Drouard</b> Complete non modular instrument</p>		<p>Angled handle</p>					
Jaws Scale 1/1		<b>Modular</b>			Short	Standard	Large
		Complete instruments					
					<p>CEV8445IU 235 mm</p>	<p>CEV445IU 335 mm</p>	<p>CEV9445IU 435 mm</p>
<p><b>Drouard</b> Modular instrument</p>		<p>Rotatable 40 -&gt; 45</p>	<p>Non rotatable 46 - 47</p>				

<b>Modular</b>				
Inserts				Handle
Blades Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Straight, serrated blades 10mm	<b>CEV8104IU</b> 240 mm	<b>CEV104IU</b> 340 mm	<b>CEV9104IU</b> 440 mm	
 Straight, Double-action jaws, 10mm	<b>CEV8623T1U</b> 240 mm	<b>CEV623T1U</b> 340 mm	<b>CEV9623T1U</b> 440 mm	
 Curved, Double-action jaws, 15mm	<b>CEV8605T1U</b> 245 mm	<b>CEV605T1U</b> 345 mm Pack 1  <b>CEV605T3U</b> 345 mm Pack 3	<b>CEV9605T1U</b> 445 mm	  
 Curved, double-action jaws, 15mm with insulated tube delivered with CEV39L (see page 71)	-	<b>CEV605T5U</b> 345 mm Pack 3	-	<b>Rotatable</b> 40 -> 45
 Curved, insulated, double-action jaws, 15mm	<b>CEV8605GU</b> 250 mm Pack 3	<b>CEV605GU</b> 350 mm Pack 3	<b>CEV9605GU</b> 450 mm Pack 3	<b>Non rotatable</b> 46 - 47
 Curved, insulated, double-action titanium jaws, 15mm	<b>CEV8605GTU</b> 250 mm	<b>CEV605GTU</b> 350 mm	<b>CEV9605GTU</b> 450 mm	  
 Curved, Double-action jaws, 10mm	<b>CEV8604T1U</b> 240 mm	<b>CEV604T1U</b> 340 mm Pack 1  <b>CEV604T3U</b> 340 mm Pack 3	<b>CEV9604T1U</b> 440 mm	






<b>Modular</b>				
Insulated insert	Wire insert			Handle
Double-action jaws for wire insert Scale 1/1	Working length - Scale 1/6			Page #
	245 mm	345 mm	445 mm	
 <b>CEV607U</b> Curved, Metzenbaum type, 15 mm - Pack 6	  			  <b>Rotatable</b> 40 -> 45
 <b>CEV607CU</b> Curved blades, 10mm - Pack 6	Wire insert without insulated sheath - single <b>CEV8607T1U</b> <b>CEV607T1U</b> <b>CEV9607T1U</b>			
 <b>CEV607DU</b> Straight blades, 10mm - Pack 6				
 Curved, Metzenbaum type <b>CEV605BGPU</b> single <b>CEV605GPU</b> Pack 3	  			<b>Non rotatable</b> 46 - 47  
à l'unité Pack 3	<b>CEV8649BGPU</b>	<b>CEV649BGPU</b>	<b>CEV9649BGPU</b>	
	<b>CEV8649GPU</b>	<b>CEV649GPU</b>	<b>CEV9649GPU</b>	
	Black sheath for metal tube - Pack 3			
	<b>CEV8649C5U</b>	<b>CEV649C5U</b>	<b>CEV9649C5U</b>	

**Nota:**  
For cleaning, the insert and the wire insert must be disassembled.



<b>NON Modular</b>				
<b>Complete instruments</b>				
Jaws Scale 1/1	Handle	Working length		
		Short	Standard	
 Straight scissors, serrated blades, 10mm	 Angled	<b>Montupet</b> <b>CEV804</b> 235 mm	<b>Manhès</b> <b>CM104</b> 355 mm	
	 Axial	<b>CEV904</b> 235 mm	–	
 Straight scissors delicate serrated blades, 10mm	 Angled	–	<b>Manhès</b> <b>CM106</b> 355 mm	
 Curved scissors, Serrated blades, 10mm	 Angled	–	<b>CEV114MU</b> 355 mm	
	 Axial	–	<b>CEV514MU</b> 335 mm	





## NON Modular





### Complete instruments

Jaws Scale 1/1	Handle	Working length
 <p>Curved scissors Serrated blades, 25 mm</p>	 <p>Angled</p>	<p><b>CEV404</b> 340 mm</p>
	 <p>Axial</p>	<p><b>CEV503</b> 350 mm</p>

## Modular


 <p>Curved scissors, serrated blades, Double-action jaws 25 mm</p>	 <p><b>Rotatable</b> 44 -&gt; 45</p>	<p><b>CEV653-1</b> 365 mm</p>
---	--	-----------------------------------




<b>Modular</b>				
Inserts			Handle	
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Straight clamp insert 25mm jaws	CEV8626I 250 mm	CEV626-1 350 mm	-	 <b>Rotatable</b> 40 -> 45
 Curved clamp insert Double-action jaws 25mm jaws	CEV8627I 255 mm	CEV627I 355 mm	CEV9627I 455 mm	<b>Non rotatable</b> 46 - 47 

<b>NON Modular</b>		
Complete instruments		
Jaws Scale 1/1	Handle	Working length
		Standard
 <b>Debakey</b> Straight clamp 20mm jaws	 Axial, with ratchet	<b>Debakey</b> CEV526 350 mm
 <b>Desmaizieres</b> Straight clamp 20mm jaws	 Angled, with long ratchet	<b>Desmaizieres</b> CEV527 365 mm

**NON Modular**

**Complete instruments**









	Jaws Scale 1/1	Handle	Working length
	 <p><b>Dulucq</b> Left curved, 60mm jaws</p>	 <p>Axial, with ratchet</p>	<p><b>Dulucq</b> <b>CEV570</b> 400 mm</p>





<b>Modular</b>		
Inserts		Handle
Jaws Scale 1/1	Working length	Page #
 <p><b>Delaitre</b> Left curved, 45mm jaws</p>	<p>CEV635-1 385 mm</p>	 <p><b>Rotatable</b> <b>44 -&gt; 45</b></p> 



**NON Modular**

**Complete instruments**

Jaws Scale 1/1	Handle	Working length
 <p><b>Querleu</b> Left curved, 45mm jaws</p>	 <p>Angled with long ratchet</p>	<p><b>Querleu</b> CEV435-1 375 mm</p>
 <p><b>Delaitre</b> Left curved, 45mm jaws</p>	 <p>Axial, with ratchet with safety lock</p>	<p><b>Querleu</b> CEV535 385 mm</p>
 <p><b>Querleu</b> Right curved, 45mm jaws</p>	 <p>Angled with long ratchet</p>	<p><b>Delaitre</b> CEV505 385 mm</p>
 <p>Curved jaws 75mm jaws</p>	 <p>Axial, with ratchet with safety lock</p>	<p><b>Querleu</b> CEV435-2 375 mm</p>

<b>Modular</b>			
Inserts		Handle	
Jaws Scale 1/1	Working length	Page #	
 <p style="text-align: center;">Straight clamp insert 65x7mm jaws Ø 10mm</p>	<p>CEV608-1 400 mm</p>		
 <p style="text-align: center;"><b>Annibale</b> Straight clamp insert Double-action jaws 65x7mm jaws Ø 10mm</p>	<p><b>Annibale</b> CEV608D 400 mm</p>	<p><b>Rotatable</b> 44 - 45</p> 	














**NON Modular**

**Complete instruments**








Jaws Scale 1/1	Handle	Working length
	 <p>Angled with long ratchet</p>	<p><b>Delaitre</b> <b>CEV431-65</b> 390 mm</p>
<p>Straight clamp, 65x7mm jaws Ø 10mm</p>	 <p>Axial, with ratchet with safety lock</p>	<p><b>CEV508</b> 400 mm</p>





NON Modular			
Complete instruments			
Scale 1/1	Jaws	Handle	Reference
	 <p><b>Satisfy</b> Angled clamp, curved tube, working length 260 mm, 110 mm jaws</p>	 <p>Axial, with ratchet with safety lock</p>	CEVA101
	 <p>Angled clamp, curved tube, working length 250 mm, 90 mm jaws</p>		CEVA102
	 <p><b>Satisfy</b> Angled clamp, curved tube, working length 240 mm, 80 mm jaws</p>		CEVA103
	 <p>Angled clamp, curved tube, working length 290 mm, 110 mm jaws</p>		CEVA104 290 mm

**NON Modular**




**Complete instruments**




Scale 1/1	Jaws	Handle	Working length
	 <p><b>Desmaizieres</b> Angled clamp, curved tube, working length 270 mm, 85 mm jaws</p>		CEVA110
	 <p><b>Desmaizieres</b> Angled clamp, curved tube, working length 270 mm, 90 mm jaws</p>	 <p>Axial, with ratchet with safety lock</p>	CEVA528
	 <p><b>Desmaizieres</b> Angled clamp, curved tube, working length 220 mm, 60 mm jaws</p>		CEVA529

Temporary Clamps

 <p>CEV565 + CEV560</p>		
<p>Jaws Scale 1/1</p>		
	<p>Straight vascular clamp 15x3mm Pack 2</p>	<p><b>Sartori</b> CEV568</p>
	<p>Straight vascular clamp 25x3mm Pack 2</p>	<p><b>Begin</b> CEV565</p>
	<p>Choledoque clamp 35x4mm Pack 2</p>	<p><b>Begin</b> CEV566</p>
	<p>Intestinal clamp 45x5mm Pack 2</p>	<p><b>Begin</b> CEV567</p>

Clamp applicator, Ø 10 mm


<b>NON Modular</b>			
Jaws Scale 1/1	Handle		Working length
  <p>Clamp applicator 360mm</p>		<p>Axial, with ratchet</p>	<p><b>Begin</b> CEV560</p>

<b>Modular</b>		
  <p>Clamp applicator insert for tube Ø10mm 310mm</p>	 <p><b>Rotatable</b> 44 -&gt; 45</p>	<p>CEV560-1 355 mm</p>

**Clips applicator**






**NON Modular**









**Complete instruments**

Jaws Scale 1/1	Handle		Working length
 Clip applicator for Tl300 et LC300 clips Ø 10mm		Angled with spring	<b>CEV310</b> 335 mm

**Circular suturing holder insert**

**Modular**

 <b>Morino</b> Circular suturing holder insert Double-action jaws Ø 10mm	 <b>Rotatable</b> 44 -> 45	<b>Morino</b> CEV603I 375 mm
 <b>Massard</b> Circular suturing holder insert Ø 5mm	 <b>Rotatable 40 -&gt; 45</b>  <b>Non rotatable 46 - 47</b>  	<b>Massard</b> CEV622IU 335 mm

<b>Modular</b>			
<b>Rotatable Instruments</b>			
Angled handles with black rotatable ring Scale 1/3			
 with ratchet	<p><b>CEV1019T5BU</b> Large pattern</p>	 with ratchet	<p><b>CEV10195CU</b></p>
 with ratchet	<p><b>CEV10195DU</b></p>	 with ratchet	<p><b>CEV10195RU</b></p>
 without ratchet	<p><b>CEV1039T5TBU</b> Large pattern</p>	 without ratchet	<p><b>CEV10395DU</b></p>
 without ratchet	<p><b>CEV10395RU</b></p>	 Insulated handle without ratchet	<p><b>CEV1039G5U</b></p>

Tube for angled handles, 5 et 3.5mm  
see page 43

## Modular

### Rotatable Instruments

Axial handles with black rotatable ring  
Scale 1/3



CEV6285R

with ratchet



CEV628P5

with ratchet



CEV6385RU





without ratchet



CEV728T5

with long ratchet

Tube for rotatable axial handles, 5 et 3.5mm  
see page 43





<b>Modular</b>	
<b>Rotatable Instruments</b>	
Ergonomic handles with black rotatable ring Scale 1/3	
 Large ergonomic handle with ratchet  <b>CEV720BRUC</b>	 Fine ergonomic handle with ratchet  <b>CEV720RU</b>
 <b>Elhadad</b> Fine ergonomic handle with ratchet  <b>CEV7285RU</b>	 Large ergonomic handle with ratchet  <b>CEV728BR5U</b>

Tube for rotatable ergonomic handles, 5 et 3.5mm  
see page 43





**Insulated tubes - 5 mm**



<b>Modular</b>		
<b>Rotatable instruments</b>		
Working length		
Short 210mm	Standard 310 mm	Long 410mm
		
Black insulated tube		
CEV8649T5TBU	CEV649T5BU	CEV9649T5TBU
		
Yellow insulated tube		
CEV8646B5U	CEV646B5U	CEV9646B5U
		
Blue insulated tube		
CEV8647B5U	CEV647B5U	CEV9647B5U
		
Green insulated tube		
CEV8648B5U	CEV648B5U	CEV9648B5U

**Metal tube and Pack of 3 black sheaths**










		
Metal tube 5mm diameter for black sheath		
CEV8649B5U	CEV649B5U	CEV9649B5U
		
Black sheath 5mm diameter for tube Pack 3		
CEV8649C5U	CEV649C5U	CEV9649C5U



## Modular

### Rotatable Instruments

Handles with grey rotatable ring  
Scale 1/3

 <p>Angled handle with long ratchet</p>	 <p>Angled handle without ratchet</p>
 <p>Axial handle with ratchet</p>	 <p>Axial handle without ratchet</p>
 <p>Fine ergonomic handle with ratchet</p>	 <p>Large ergonomic handle with ratchet</p>
 <p>Axial handle with long ratchet</p>	<p>Tube for rotatable handles for <math>\varnothing</math> 5 and 10 mm non-insulated tube see page 45</p>



## Modular

### Rotatable instruments

Working length

Short

Standard

Large



Non-insulated tube for rotatable instrument  $\varnothing$  10mm

CEV864910B  
215mm

CEV649-10B  
305mm

–








Non-insulated tube for rotatable instrument  $\varnothing$  5mm

CEV8649X5  
210mm

CEV649X5  
310mm

CEV9649X5  
410mm

**Non Rotatable handles - 5 mm - for black insulated tube**

<b>Modular</b>	
<b>Non rotatable Instruments</b>	
Handles without rotatable ring Scale 1/3	
 Angled handle with long ratchet CEV1019T5NUC	 Axial handle with ratchet CEV628-5N
 Angled handle with long ratchet CEV10195NAU	 Axial handle without ratchet CEV638T5NU
 Angled handle without ratchet CEV1039T5NU	

**Black insulated tube for non-rotatable instrument  $\varnothing$  5mm**

Black insulated tube for non-rotatable instrument $\varnothing$ 5mm		
Working length		
Short 210mm CEV8649T5NU	Standard 310mm CEV649T5NU	Large 410mm CEV9649T5NU

**Non Rotatable handles - 5 mm - for non-insulated tube**

**Modular**

**Non rotatable instruments**


**Handles without rotatable ring**  
Scale 1/3



 <p>CEV720</p> <p>Fine ergonomic handle with on/off ratchet</p>	 <p>CEV720B</p> <p>Large ergonomic handle with on/off ratchet</p>
--	---

**Non-insulated tube for non-rotatable instrument ø 5mm**



 <p>Non-insulated tube for non-rotatable instrument ø 5mm</p>		
<p>Working length</p>		
<p>Short 210mm</p>	<p>Standard 310mm</p>	<p>Large 410mm</p>
<p>CEV8649P5</p>	<p>CEV649P5</p>	<p>CEV9649P5</p>

Bipolar cables 4.5 m



**CP393-2**

Bowa Emc Lamidey Valleylab













**CP393-3**

Erbe VIO ICC ACC  
(except international)  
Storz Wolf











**CP393-4**

Aesculap Berchtold Martin  
Sutter BM 780 II









<b>Modular</b>								Handles Non rotatable
Inserts & Tubes								
Jaws Scale 1/1	Working length							
	Short		Standard		Large			
	Insert	Tube 1/3	Insert	Tube 1/3	Insert	Tube 1/3		
	 <b>Mouiel</b> 15x5mm	CEV86341A 245 mm		CEV634-1A 340 mm		CEV9634-1A 430 mm	 CEV669B Angled handle	
	 with hook 15x5mm	-		CEV632-1A 340 mm		CEV96321A 430 mm	 CEV669BG <b>Gayet</b> Ergonomic handle	
	 <b>Gayet</b> 20x5mm	CEV86361A 250 mm		CEV636-1A 345 mm		CEV96361A 435 mm	 MCEN54P <b>Dessi</b> Angled handle	
	 <b>Vallancien</b> with teeth 20x5mm	-		CEV6371 345 mm		-	 CEV86795B 210 mm	 CEV6795B 310 mm












<b>NON Modular</b>					
Complete instruments					
Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
 15x5mm	 Angled	CEV834 240 mm	<b>Mouiel</b> CEV134 360 mm	–	
	 Ergonomic with graduated tube 100mm	–	<b>Gayet</b> CEV134BG 360 mm	<b>Gayet</b> CEV9134BG 450 mm	
	 Axial	–	CEV580-2 365 mm	–	
 15x5mm with hook	 Angled	–	CEV132 360 mm	–	
 <b>Gayet</b> 20x5mm	 Angled	–	<b>Gayet</b> CEV136 365 mm	–	





<b>Modular</b>					Handles Non rotatable
Inserts & Tubes					
Jaws Scale 1/1	Working length				
	Standard		Large		
	Mors	Tube 1/3	Mors	Tube 1/3	
 <b>Manhès</b> Spatula jaws 10x3mm	CEV630-1A 340 mm		CEV9630-1A 425 mm		 <b>CEV669B</b> Angled handle
 <b>Botella</b> Triangle fenestrated jaws 12x5mm	CEV633-1A 340 mm		CEV96331A 430 mm		 <b>Gayet</b> CEV669BG Ergonomic handle
 Spatula angled jaws 10x2mm	CEV630C 340 mm		-		 <b>MCEN54P</b> <b>Dessi</b> Angled handle
		CEV6795B 310 mm		CEV96795B 395 mm	



<b>NON Modular</b>				
Complete instruments				
Jaws Scale 1/2	Handle		Working length	
			Short	Standard
 Spatula jaws 10x3mm		Angled	<b>Montupet</b> CEV830 240 mm	<b>Manhès</b> CEV130 360 mm
 <b>Botella</b> Triangle fenestrated jaws 12x5mm		Angled	—	<b>Botella</b> CEV133 360 mm
		Axial	—	<b>Botella</b> CEV580-3 360 mm
—	—	—	—	—
 Delicate triangular jaws 8x5mm		Angled	—	<b>Botella</b> CEV133-1 360 mm
 <b>Manhès</b> 10x2mm		Angled	—	<b>Manhès</b> CM131 360 mm



## Modular

### Hook jaws

Scale 1/1



**CEV229-2A**  
Crochet # 2  
Pack 3



**CEV229-3A**  
Crochet # 3  
Pack 3



**CEV229-4A**  
Crochet # 4  
Pack 3



**CEV2295A**  
Crochet # 5  
Pack 3

### Manche de crochet

Scale 1/3

**CEV225A**

Working length 430mm



**CEV82291**

Working length 230mm



**CEV229-1A**

Working length 350mm



**CEV9229-1**

Working length 450mm



**CEV229B**

Working length 380mm



Page 54

## Monopolar hook, non modular

Scale 1/3



1/1



**CEV225**

Working length 430mm



1/1



**Périssat**  
**CEV230-1**

Working length 330mm







1/1





**Montupet**  
**CEV825**

Working length 220mm

**Retractable knife, monopolar electrode**

		<b>NON Modular</b>	
		Complete instruments	
1/3			
1/1		<b>CEV135-5</b>	Retractable knife ø 5mm 300mm without knife blades SM63 / SM64 / SM65
1/3			
1/1		<b>CM109</b> <b>Manhès</b>	Retractable monopolar electrode ø 5mm 320mm









**Monopolar palpation probe**

1/3			
1/1		<b>CEV443</b> <b>Delaitre</b>	Monopolar palpation probe ø 5mm 350mm














<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex

Palpation probe

	1/3				
	1/1		<b>CEV204</b> 320 mm - Ø 5mm	Round probe	
	1/1		<b>CEV204G45</b> 450 mm - Ø 5mm	Round probe, graduated 10 to 150 mm	
	1/1		<b>CEV205</b> 340 mm - Ø 5mm	Flat probe	
	1/3				
	1/1		<b>CEV207</b> 410 mm - Ø 5mm	Double ended probe, round end flap graduated 10 to 150 mm	

Articulated retractor

			<b>Flexible retractor</b>  <b>CEV416</b> 380mm - Ø 5mm		1/3
			<b>Retractor with 3 blades 35mm</b>  <b>CEV413</b> 370mm - Ø 5mm		1/3
	<b>Nouaille</b>		<b>Retractor diamond design</b>  <b>CEV425</b> 320mm - Ø 5mm		1/3
			<b>Retractor with 3 blades 65mm</b>  <b>CEV423</b> 400mm - Ø 10mm		1/3
			<b>Retractor with 5 blades 65mm</b>  <b>CEV453</b> 400mm - Ø 10mm		1/3
			<b>Articulated retractor with 5 blades 65mm</b>  <b>CEV453A</b> 330mm - Ø 10mm		1/3

<b>NON Modular</b>	
<b>Complete instruments</b>	
Flexible tension arm delivered with instruments holder CEV980-3	Accessories for flexible tension arm
 <p><b>Mosnier</b> 500mm : CEV980A 600mm : CEV980B delivered with instruments holder CEV980-3</p>	
	<p>CEV980-3 Instruments holder Ø 5 et 10mm <b>Mosnier</b></p> 
	<p>CEV984 Operating table clamp with screw</p> 
	<p>CEV985 Operating table clamp</p> 



### CM201 irrigation suction pump accessories

Suction irrigation tubes  
Sterile - Pack 10

Electric handle

CM217S<sup>(D)</sup>

CM201-9



### Eyes suction irrigation tubes for handle CM201-9

1/1

«L» Working length

1/3

Eyes hydro cannula Ø 10mm



CM20132H10  
L 320mm  
CM20128H-10  
L 280mm



Eyes hydro cannula Ø 5mm



CM20132H-5  
L 320mm  
CM20128H-5  
L 280mm  
CM20140H5  
L 400mm



Eyes cannula Ø 5mm



CM20132L-5  
L 320mm



### Cleaning sleeve for telescope




Cleaning sleeve for telescope odeg  
Ø 10mm 300mm  
Scale 1/3

CEV81

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex

Suction irrigation tube








Modular suction irrigation tube


Eyes suction irrigation tube Ø 5mm		Handles for tube Ø 5mm
<p>1/1</p> <p>195mm CEV8103D5</p> <p>295mm CEV103-5</p> <p>395mm CEV9103D5</p>	<p>1/3</p> <p>CEV103AH</p>	
<p>1/1</p> <p>Blunt tip</p> <p>195mm CEV8103G5</p> <p>295mm CEV103G5</p> <p>395mm CEV9103G5</p>	<p>1/3</p> <p>CEV103BG</p> <p>1/3</p> <p>CEV103GH</p>	
Eyes suction irrigation tube Ø 10mm		Handles for tube Ø 10mm
<p>1/1</p> <p>CEV103-10</p> <p>295mm</p> <p>1/3</p>	<p>1/3</p> <p>CEV103AH</p>	
<p>1/1</p> <p>CEV103P</p> <p><b>Prete</b></p> <p>410mm</p>		









**Wire-pass**

	<p><b>Articulated</b></p> <p><b>CEV438</b></p> <p>420mm - Ø 10mm</p>	<p><b>Bertrand</b></p>		<p>1/3</p>
	<p><b>Articulated</b></p> <p><b>CEV438B</b></p> <p>335mm - Ø 5mm</p>			<p>1/3</p>
	<p><b>Retractable</b></p> <p><b>CEV438D5</b></p> <p>430mm - Ø 5mm</p>	<p><b>Morino</b></p>		<p>1/3</p>
	<p><b>Retractable</b></p> <p><b>CEV456</b></p> <p>170mm - Ø 2.5mm</p>			<p>1/3</p>

**Mintz set for punction**


<p>Sleeve - Ø 2.5mm - 240mm</p> <p><b>CM170-1</b></p>		<p>1/3</p>
<p>Sharp trocar</p> <p><b>CM170-2</b></p>		
<p>Smooth trocar</p> <p><b>CM170-3</b></p>		

**Needles**

<p>Insufflation needle with luer-lock</p> <p>120mm</p> <p><b>CEV60-120</b></p>		<p>1/3</p>
<p>150mm</p> <p><b>CEV60-150</b></p>		
<p>Punction needle with luer-lock</p> <p>Ø 5mm - 300mm</p> <p><b>CEV65</b></p>		
<p>Ø 1.2mm - 40mm</p> <p><b>CEV24F</b></p>		
<p>Ø 2.0mm - 40mm</p> <p><b>CEV24G</b></p>		<p>with a lateral perforation to purge</p> <p> 2/1</p>

**Cholangiography**


1/3



1/1



1/3




Sealing cap

**CEV440P6**

Pack 6



1/3



**Falchero**



Sealing cap

**CEV103P6**

Pack 6

**Cholédoscopy**

1/3



1/3












Sealing cap


**CEV105P6**

Pack 6



**Trocar with metal cannula**








			Working length	
			50 mm	80 mm
Trocar tip	smooth		CEV8143A4	CEV143A4
	pyramidal		CEV8143A3	CEV143A3
Sealing cap	green		CEV103P6 Pack 6	
Trocar head	magnetic		CEV177A3 Delivered complete with 1 magnetic valve and 1 green sealing cap	
	silicone		CEV127T3 Delivered complete with 1 silicone valve and 1 green sealing cap	
Valve for trocar head	magnetic		CEV179A5 Pack 6	
	silicone		CEV125P6 Pack 6	
Metal cannula, crenelated	with stopcock		CEV8143A7	CEV143A7
	without stopcock		CEV8143B7	CEV143B7



Scale 1/1

Composite trocar with silicone valve








			Working length	
			45 mm	70 mm
Trocar tip 	smooth 		CEV8826M	CEV326M
	pyramidal 		CEV8826P	CEV326P
Sealing cap	green 		CEV103P6 Pack 6	
Valve for trocar head	silicone 		CEV125P6 Pack 6	
Trocar body composed of 1 trocar head with composite cannula and silicone valve	smooth cannula without stopcock 		CEV8826L	CEV326L
	threaded cannula without stopcock 		CEV8826V	CEV326V





Scale 1/1



CEV8126V

			Working length			
			50 mm	80 mm	100 mm	120 mm
Trocar tip 	blunt 		CEV8145A5	CEV145A5	CEV7145A5	CEV9145A5
	sharp 		CEV8145A3	CEV145A3	CEV7145A3	CEV9145A3
Sealing cap	blue 		CEV105P6 Pack 6			
Head of trocar	magnetic valve 		CEV177A5 delivered with one magnetic valve and one blue sealing cap			
	silicone valve 		CEV127T delivered with one silicone valve and one blue sealing cap			
Valve for head of trocar	magnetic 		CEV179A5 Pack 6			
	silicone 		CEV125P6 Pack 6			
Sleeve	with stopcock	screw 	CEV8145A6	CEV145A6	CEV7145A6	CEV9145A6
		smooth 	CEV8145A7	CEV145A7	CEV7145A7	CEV9145A7
	without stopcock	screw 	CEV8145B6	CEV145B6	CEV7145B6	CEV9145B6
		smooth 	CEV8145B7	CEV145B7	CEV7145B7	CEV9145B7

Body of trocar 5 mm diameter with silicone valve with screw sleeve, without stopcock 	1/2 CEV126V	CEV8126V Working length 45 mm <b>Montupet</b>	CEV126V Working length 75 mm
Trocar tip for body of trocar 5mm diameter 	blunt 	CEV8126PM Working length 95 mm <b>Montupet</b>	CEV126PM Working length 120 mm
	sharp 	CEV8126P Working length 95 mm	CEV126P Working length 120 mm



Scale 1/1

				Working length	
				90 mm	130 mm
Trocar tip		blunt		CEV150-4	CEV9150A4
		sharp		CEV150-3	CEV9150A3
Sealing cap		red		CEV110P6 Pack 6	
Head of trocar		magnetic valve		CEV177A10 delivered with one magnetic valve and one red sealing cap	
		silicone valve		CEV127T10 Delivered with one silicone valve and one red sealing cap	
Valve for head of trocar		magnetic		CEV179A10 Pack 6	
		silicone		CEV125Q6 Pack 6	
Sleeve	with stopcock	screw		CEV150A6	CEV9150A6
		smooth		CEV150A7	CEV9150A7
	without stopcock	screw		CEV150B6	CEV9150B6
		smooth		CEV150B7	CEV9150B7

### Dulucq smooth trocar tip palpator probe

Scale 2/3

Dulucq















CEV480








Ø 10 mm

Working length 180mm



Scale 1/1

			Working length	
			90 mm	130 mm
Trocar tip 	blunt 		CEV172A4	CEV9172A4
	sharp 		CEV172A3	CEV9172A3
Sealing cap	Orange 		CEV112P6 Pack 6	
Head of trocar	magnetic valve 		CEV177A12 delivered with one magnetic valve and one orange sealing cap	
	silicone valve 		CEV127T12 delivered with one silicone valve and one orange sealing cap	
Valve for head of trocar	magnetic 		CEV179A12 Pack 6	
	silicone 		CEV125Q6 Pack 6	
Sleeve	with stopcock	screw 	CEV172A6	CEV9172A6
		smooth 	CEV172A7	CEV9172A7
	without stopcock	screw 	CEV172B6	CEV9172B6
		smooth 	CEV172B7	CEV9172B7

	Ø 5.5 mm			Ø 10.5 mm		Ø 12.5 mm	
1/3	CEV805I	CEV105	CEV105I	CEV110I	CEV112I	CEV122I	CEV212
							
	Working length 50mm with stopcock	Working length 80mm without stopcock	Working length 80mm with stopcock	Working length 80mm with stopcock	Working length 120mm with stopcock	Working length 80mm with stopcock	Working length 120mm without stopcock
	<b>Montupet</b>						


### Spring for trumpet valve trocar



CEV105R  
Ø 5.5mm



CEV110R  
Ø 10.5mm



CEV112R  
Ø 12mm

Sealing cap for trumpet valve trocars  
see page 68





			Diameter		
			Ø 5.5mm	Ø 10mm	Ø 12.5mm
Trocar tip		blunt 	CEV145-4	CEV150-4	—
		sharp 	CEV145-3	CEV150-3	CEV172-3
Sealing cap			 CEV105P6 Pack 6	 CEV110P6 Pack 6	 CEV112P6 Pack 6
Head of trocar	with stopcock		 CEV178-5	 CEV278-10	 CEV278-12
	without stopcock		 CEV177-5	—	 CEV177-12
Sleeve	without stopcock	screw	 CEV145-6	 CEV150-6	 CEV172-1
		smooth	 CEV145-7	 CEV150-7	—

Sealing cap for trocar



CEV103P6

Ø 3.5mm

Pack 6



CEV105P6

Ø 5.5mm

Pack 6



CEV110P6

Ø 10.5mm

Pack 6



CEV112P6

Ø 12mm

Pack 6

Reducer sleeve for trocar

1/1



CEV155A

from Ø 10.5 to 5.5mm



CEV156A

from Ø 12/13 to 5.5mm



CEV181A

from Ø 12/13 to 10.5mm

Silicone ring for attaching 15mm reducer sleeve



CEV155B

Pack 6

1/3



CEV153

135 mm

from Ø 5.5 to 3.5mm



CEV155

135 mm

from Ø 10.5 to 5.5mm



CEV155-1

65 mm

from Ø 10.5 to 5.5mm



CEV155P

135 mm

from Ø 12/13 to 5.5mm

Tube interne Ø 5mm



CEV156

135 mm

from Ø 12/13 to 10.5mm



CEV181

135 mm

from Ø 12/13 to 10.5mm



CEV181-1

65 mm

from Ø 12/13 to 10.5mm

Double retractor

1/3



**Descottes**

CEV417

200 mm

pair

Double retractor, blunt



CEV418

200 mm

pair

Double retractor, blunt



**Mosnier**










CEV419

160 mm








pair

Double retractor, blunt




**Trocars 85mm**

Scale 1/3	Blunt trocar tip	Sharp trocar tip	Flexible threaded cannula
	<b>CET105-1</b>	<b>CET105-2</b>	<b>CET105-3</b>
Ø 5.5mm			
	<b>CET110-1</b>	<b>CET110-2</b>	<b>CET110-3</b>
Ø 10.5mm			
	<b>CET112-1</b>	<b>CET112T2</b>	<b>CET112-3</b>
Ø 12.5mm			
Joints : See page 68			

**Forceps, non modular instruments**

Jaws Scale 1/1	Handle	Reference
 <b>Manhès</b> curette forceps 12x5mm Working length 355mm	 Angled	<b>CM107U</b>
 curette forceps 12x5mm, Double-action jaws Working length 355mm	 Angled	<b>CEV444</b>
 fenestrated forceps 20x5mm, curved tube Working length 340mm	 Axial, with ratchet   Curved tube	<b>CET525</b>


**Modular forceps**

	 <b>Manhès</b> curette insert 12x5mm Ø 5mm	Working length		
		Short	standard	Large
 <b>Rotable : page 40 - 45</b>   <b>Non rotatable : page 46 - 47</b>	<b>CEV8107IU</b> 245mm	<b>CEV107IU</b> 345mm	<b>CEV9107IU</b> 445mm	



Belghiti	Soubrane	CEV392P20	Dousset	Gayet	Partenski
215mm	215mm	215mm	290mm	220mm	190mm
CEV392D20	CEV392S20	CEV392P20	CEV392D30	CEV392G24	CEV392R15
straight	straight	straight	straight	straight	bayonet
Jaws 2mm	Jaws 2mm double irrigation	Serrated jaws 2mm	Jaws 3mm double irrigation	Jaws 5x15mm	Jaws 4x10mm

Bipolar cables - 4,5 m



**CP393-2**  
Bowa Emc Lamidey Valleylab














**CP393-3**  
Erbe VIO ICC ACC (except international) Storz Wolf




**CP393-4**  
Aesculap Berchtold Martin Sutter BM-78o II












**Semi-modular**

Insert with tubes		Handle
Insert insulated tube Scale 1/1	Working length	Page #
 Babcock Insert	CEV391IU 335 mm	      <b>Rotatable</b> <b>40 -&gt; 45</b>
 Curved scissors insert, serrated blades 8 mm	CEV393IU 330 mm	
 Straight scissors insert, serrated blades 8mm	CEV394IU 330 mm	
 Fenestrated insert, 3x12mm jaws	CEV395IU 335 mm	
 Alligator jaws insert 6 mm	CEV396IU 330 mm	
 Dissector insert, medium curved 10 mm jaws	CEV398IU 340 mm	
 <b>Allis</b> Grasping insert	CEV399IU 335 mm	
-	-	
-	-	



CEV39L  
Adapter cleaning luer-lock for tube insert












<b>NON Modular</b>			
<b>Complete instruments</b>			
Jaws Scale 1/1	Handle		Working length
 Babcock forceps	 Angled with long ratchet		<b>CEV391BU</b> 355 mm
-	-	-	-
 Straight scissors, serrated blades, 8mm	 Angled		<b>CEV394BU</b> 350 mm
 Fenestrated forceps, 3x12mm jaws	 Angled with long ratchet		<b>CEV395BU</b> 355 mm
-	-	-	-
-	-	-	-
 <b>Allis</b> Grasping forceps with teeth	 Angled with long ratchet		<b>CEV399BU</b> 355 mm
 Needle-holder, tungsten straight jaws	 Axial		<b>CEV992B</b> 340 mm
 Needle-holder, tungsten left curved jaws			<b>CEV992C</b> 340 mm

**NON Modular**

**Complete instruments : Bipolar forceps**



Jaws Scale 1/1	Handle	Working length
 <p><b>Gayet</b> Forceps, fenestrated jaws, 10x3mm, graduated tube 100mm</p>	 <p>Ergonomic bipolar</p>  <p>1/3</p>	<p><b>CEV397BG</b> 305mm</p>
 <p>Forceps, 2mm spatula jaws</p>	 <p>Ergonomic bipolar</p>  <p>1/3</p>	<p><b>CEV397C</b> 305mm</p>
 <p>Forceps, 2mm serrated jaws</p>	 <p>Angled bipolar</p>  <p>1/3</p>	<p><b>CEV397</b> 305mm</p>

**Cadière forceps**



Forceps without port  
Fenestrated forceps  
percutaneous


CEV350

1/1

195mm



Cadière



1/3

**Palpation probe**



Round probe Ø 3.5mm 330mm

CEV203

1/1



1/3

**Monopolar hooks**



1/1



CEV390B

Ø3.5mm, Working length 330mm, round tip



1/3

1/1



CEV390C

Ø3.5mm Working length 355mm, 90°



1/3



Page 54

**Suction irrigation tube**



Eyes suction irrigation tube, graduated 50/100/150mm

1/1



CEV103G3

Ø 3.5mm 295mm



1/3

Graduated suction irrigation tube for handle CEV103GH or CEV103BG

Handle for suction tube CEV103G3



CEV103GH



CEV103BG

1/3

**Iannelli knot pusher**



1/1



CEV342

440mm working length 340mm Ø 3.5mm



1/3

1/1



CEV342E









360mm working length 340mm Ø 3.5mm




1/3



**Trocart métal**

			Working length	
			50 mm	80 mm
Trocar tip	blunt		CEV8143A4	CEV143A4
	sharp		CEV8143A3	CEV143A3
Sealing cap	Green		CEV103P6 Pack 6	
Head of trocar	magnetic valve		CEV177A3 delivered with one magnetic valve and one green sealing cap	
	silicone valve		CEV127T3 delivered with one silicone valve and one green sealing cap	
Valve for head of trocar	magnetic		CEV179A5 Pack 6	
	silicone		CEV125P6 Pack 6	
Metal sleeve	with stopcock		CEV8143A7	CEV143A7
	without stopcock		CEV8143B7	CEV143B7










Scale 1/1








Trocart with composite sleeve and silicone valve








Scale 1/1

			Working length	
			45 mm	70 mm
Trocar tip 	blunt 		CEV8826M	CEV326M
	sharp 		CEV8826P	CEV326P
Sealing cap	green 		CEV103P6 Pack 6	
Valve for head of trocar	silicone 		CEV125P6 Pack 6	
Body of trocar with head of trocar with silicone valve and composite sleeve	with smooth sleeve without stopcock 		CEV8826L	CEV326L
	with screw sleeve without stopcock 		CEV8826V	CEV326V








Monopolar hooks

	 <p>1/1</p>	 <p>1/3</p>
	<p><b>CEV890</b> 225mm Ø 3.5mm</p>	<p><b>Montupet</b></p>
	 <p>1/1</p>	 <p>1/3</p>
	<p><b>CEV890C</b> 90deg 225mm Ø 3.5mm</p>	 <p>Page 54</p>

Retractable knife

	 <p>1/1</p>	 <p>1/3</p>
	<p><b>CEV898</b> 90mm Ø 3.5mm delivered with one blade CEV898L</p>	
 <p>1/1</p>		
	<p><b>CEV898L</b> Knife blade for CEV898 Pack 3</p>	










Suction irrigation tube

	 <p>1/1</p>	<p>Eyes suction irrigation tube, graduated 50/100/150mm</p>  <p>1/3</p>
	<p><b>CEV8103G3</b> Ø 3.5mm 195mm</p>	<p>Eyes suction irrigation tube, blunt graduated 50/100/150mm for handle CEV103GH or CEV103BG</p>
		 <p><b>CEV103GH</b> Handle</p>
		 <p><b>CEV103BG</b> Handle</p>
 <p>1/1</p>	 <p>1/3</p>	
<p><b>CEV896</b> Ø 3.5mm 220mm</p>	<p>Eyes suction irrigation tube with blunt tip</p>	



**NON Modular**













**Complete instruments : bipolar forceps**









Jaws Scale 1/1	Handle	Working length
 <p>Bipolar forceps 1.5mm serrated jaws</p>	 <p>Angled Bipolar</p>	<p>CEV897 250 mm</p>
		
1/3		
 <p>Bipolar forceps, fenestrated jaws 10x3mm jaws</p>	 <p>Ergonomic bipolar</p>	<p>CEV897BG 250 mm</p>
		
1/3		
 <p>Bipolar forceps 2mm serrated jaws</p>	 <p>Ergonomic bipolar</p>	<p>CEV897C 250 mm</p>
		
1/3		



**Semi-Modular**



Insert with tube		Handle
Insert insulated tube Scale 1/1	Working length	Page #
 Babcock insert	CEV8911U 235 mm	      <b>Rotatable 40 -&gt; 45</b>
 Curved scissors insert, serrated blades 8 mm	CEV8931U 230 mm	
 Straight scissors insert, serrated blades 8mm	CEV8941U 230 mm	
 Fenestrated insert, 3x12mm jaws	CEV8951U 235 mm	
 <b>Allis</b> Grasping forceps	CEV8991U 220 mm	
 Alligator jaws insert, 6 mm	CEV8961U 235 mm	
 Dissector insert, Medium curved 10 mm jaws	CEV8981U 230 mm	
	CEV39L Adapter cleaning luer-lock for tube insert	









<b>NON Modular</b>			
<b>Complete instruments</b>			
Jaws Scale 1/1		Handle	Working length
 <b>Montupet</b> Babcock forceps		 Angled with long ratchet	<b>CEV891MU</b> 200 mm
—		—	—
 <b>Montupet</b> Straight scissors, serrated blades 8mm		 Angled	<b>CEV894MU</b> 220 mm
 <b>Montupet</b> Fenestrated forceps, 3x12mm jaws		 Angled with long ratchet	<b>CEV895MU</b> 195 mm
 <b>Allis</b> Grasping forceps with teeth		 Angled with long ratchet	<b>CEV899U</b> 220 mm











## NON Modular Complete instruments


### De Lagausie instruments

Jaws Scale 1/1	Handle	Working length
 <b>De Lagausie</b> Fenestrated forceps, 3x12mm jaws	 Axial	<b>CEV995MU</b> 220 mm
 <b>De Lagausie</b> Straight scissors, serrated blades 8mm	 Axial	<b>CEV996MU</b> 220 mm
 <b>De Lagausie</b> Alligator forceps, 6mm serrated jaws	 Axial	<b>CEV997MU</b> 210 mm
 <b>De Lagausie</b> Dissector, medium curved 10mm jaws	 Axial	<b>CEV998MU</b> 220 mm


### Needle-holder

 <b>Montupet</b> Needle-holder, tungsten straight jaws	 Angled with long ratchet	<b>CEV892</b> 220 mm
 <b>Reinberg</b> Needle-holder, straight jaws	 Axial	<b>CEV992</b> 220 mm
 Needle-holder, left curved tungsten jaws	 Axial	<b>CEV992D</b> 220 mm


Scale 1/3  
**Curved scissors, titanium large curved blades, 15 mm**



**Adson-Brown forceps with teeth**



**Needle-holder, tungsten straight jaws, 1 x 10 mm**



**Needle-holder, titanium automatic jaws**






Scale 1/2

Working length «L»



Trocar holder ring  
for trocar Ø 10mm


**MCLP70**



Trocar without valve 55mm  
with blunt trocar tip


**MCLP10**

L 55mm




Fenestrated forceps 3x8mm jaws

**MCLP20U**



L 180mm

Angled handle



Fenestrated forceps 3x8mm jaws

**MCLP25U**



L 100mm

Axial handle

**MCLP30U**



Forceps, alligator jaws, 6mm

**MCLP30U**



Curved scissors, blades 4mm

**MCLP40U**



L 100mm

Axial handle

**MCLP40U**



Monopolar hook

**MCLP50**



L 100mm

**MCLP50**

Page 54

Round palpation probe

**MCLP60**



L 180mm

**MCLP60**

Double blunt retractor  
180mm

**MCLT10**

Small size 8x35 / 15x40mm

**MCLT11**

large size 8x35 / 15x50mm



Blunt curved elevator  
190mm

**MCLT20**

2mm

**MCLT21**

4mm



Peanut-holder



200mm - Ø 3mm

**MCLT30**

Axial handle



Suction tube elevator

200mm - Ø 3mm

**MCLT40**

malleable



Nerve eye suction tube  
with finger control

200mm - Ø 3mm

**MCLT50**

malleable



Scheye biopsy rectal forceps, 5mm  
Hirschsprung's disease

1/3




Curette jaws insert and tube Ø 5mm, non sterile, single use (pack 3) for biopsy



Handle non sterile, resterilizable

1/4



Biopsy diameter  
3 or 5mm

CEV504S3 Insert Ø 3mm


CEV504S5 Insert Ø 5mm

CEV504PS Axial handle



**Precaution:**

If the single use symbol is on the device label then this device is designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.



Stainless steel sterilization cases

<b>7874MF-34</b>	340x200x80mm
<b>7874MF-40</b>	400x200x90mm
<b>7874MF-61</b>	610x180x80mm
<b>CEV100-3</b>	3 perforated trays for 7874MF-61



**MCO731<sup>(D)</sup>**  
Surgical instrument oil



**CEV56**  
Stainless steel cap  
for luer-lock cleaning port  
Pack 1



**CEV57**  
Silicone cap  
for luer-lock cleaning port  
Pack 6

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex

**Numerical index**

**7**

7874MF-34 .....	86
7874MF-40 .....	86
7874MF-61 .....	6

**CEC**

CEC8500T5 .....	82
CEC8522 .....	82
CEC8605GT .....	82
CEC8650T .....	82

**CET**

CET105-1 .....	69
CET105-2 .....	69
CET105-3 .....	69
CET110-1 .....	69
CET110-2 .....	69
CET110-3 .....	69
CET112-1 .....	69
CET112-3 .....	69
CET112T2 .....	69
CET525 .....	69

**CEV**

CEV24F .....	59
CEV24G .....	59
CEV39L .....	26, 71, 79
CEV56 .....	86
CEV57 .....	86
CEV60-120 .....	59
CEV60-150 .....	59
CEV65 .....	59
CEV81 .....	57
CEV100-3 .....	86
CEV101 .....	58
CEV103-5 .....	58
CEV103-10 .....	58
CEV103AH .....	58
CEV103BG .....	58, 74, 77
CEV103G3 .....	74
CEV103G5 .....	58
CEV103GH .....	58, 74, 77
CEV103P .....	59
CEV103P6 .....	60, 61, 62, 68, 75, 76
CEV104IU .....	26
CEV105 .....	66
CEV105I .....	66
CEV105P6 .....	63, 67, 68
CEV105R .....	66
CEV107IU .....	69
CEV110I .....	66
CEV110P6 .....	64, 67, 68
CEV110R .....	66
CEV112I .....	66
CEV112P6 .....	65, 67, 68
CEV112R .....	66
CEV114MU .....	28

CEV122I .....	66
CEV125P6 .....	61, 62, 63, 75, 76
CEV125Q6 .....	64, 65
CEV126P .....	63
CEV126PM .....	63
CEV126V .....	63
CEV127T .....	63
CEV127T3 .....	61, 75
CEV127T10 .....	64
CEV127T12 .....	65
CEV130 .....	52
CEV132 .....	50
CEV133 .....	52
CEV133-1 .....	52
CEV134 .....	50
CEV134BG .....	50
CEV135-5 .....	54
CEV136 .....	50
CEV143A3 .....	61, 75
CEV143A4 .....	61, 75
CEV143A7 .....	61, 75
CEV143B7 .....	61, 75
CEV145-3 .....	67
CEV145-4 .....	67
CEV145-6 .....	67
CEV145-7 .....	67
CEV145A3 .....	63
CEV145A5 .....	63
CEV145A6 .....	63
CEV145A7 .....	63
CEV145B6 .....	63
CEV145B7 .....	63
CEV150-3 .....	64, 67
CEV150-4 .....	64, 67
CEV150-6 .....	67
CEV150-7 .....	67
CEV150A6 .....	64
CEV150A7 .....	64
CEV150B6 .....	64
CEV150B7 .....	64
CEV153 .....	68
CEV155 .....	68
CEV155-1 .....	68
CEV155A .....	68
CEV155B .....	68
CEV155P .....	68
CEV156 .....	68
CEV156A .....	68
CEV172-1 .....	67
CEV172-3 .....	67
CEV172A3 .....	65
CEV172A4 .....	65
CEV172A6 .....	65
CEV172A7 .....	65
CEV172B6 .....	65
CEV172B7 .....	65
CEV177-5 .....	67
CEV177-12 .....	67
CEV177A3 .....	61, 75
CEV177A5 .....	63
CEV177A10 .....	64
CEV177A12 .....	65
CEV178-5 .....	67
CEV179A5 .....	61, 63, 75
CEV179A10 .....	64
CEV179A12 .....	65
CEV181 .....	68
CEV181-1 .....	68
CEV181A .....	68
CEV203 .....	74
CEV204 .....	55
CEV204G45 .....	55
CEV205 .....	55
CEV207 .....	55
CEV211U .....	15
CEV212 .....	66
CEV220U .....	13
CEV225 .....	53
CEV225A .....	53
CEV229-1A .....	53
CEV229-2A .....	53
CEV229-3A .....	53
CEV229-4A .....	53
CEV229B .....	53
CEV230-1 .....	53
CEV278-10 .....	67
CEV278-12 .....	67
CEV310 .....	39
CEV321 .....	15
CEV326L .....	62, 76
CEV326M .....	62, 76
CEV326P .....	62, 76
CEV326V .....	62, 6
CEV342 .....	74
CEV342E .....	74
CEV350 .....	74
CEV390B .....	74
CEV390C .....	74
CEV391BU .....	72
CEV391IU .....	71
CEV392D20 .....	70
CEV392D30 .....	70
CEV392G24 .....	70
CEV392P20 .....	70
CEV392R15 .....	70
CEV392S20 .....	70
CEV393IU .....	71
CEV394BU .....	72
CEV394IU .....	71
CEV395BU .....	72
CEV395IU .....	71
CEV396IU .....	71
CEV397 .....	73
CEV397BG .....	73
CEV397C .....	73
CEV398IU .....	71
CEV399BU .....	72
CEV399IU .....	71
CEV404 .....	29
CEV405BG .....	3
CEV405U .....	3
CEV406IU .....	4
CEV406U .....	5
CEV407IU .....	4

CEV407U.....	5	CEV525.....	3	CEV625DU.....	2
CEV410.....	9	CEV525C.....	3	CEV625EU.....	6
CEV413.....	55	CEV526.....	30	CEV625P1U.....	2
CEV416.....	55	CEV527.....	30	CEV625PU.....	6
CEV417.....	68	CEV530.....	9	CEV625T1U.....	2
CEV418.....	68	CEV531-2.....	19	CEV626-1.....	30
CEV419.....	68	CEV531T3U.....	19	CEV627I.....	30
CEV420.....	11	CEV532.....	19	CEV628-5N.....	46
CEV421.....	11	CEV532IU.....	18	CEV628-10.....	44
CEV423.....	55	CEV533-2.....	9	CEV628P5.....	41
CEV425.....	55	CEV535.....	33	CEV630-1A.....	51
CEV430.....	21	CEV538-10.....	24	CEV630C.....	51
CEV431-65.....	35	CEV538T5.....	23	CEV632-1A.....	49
CEV432.....	9	CEV538TU.....	22	CEV633-1A.....	51
CEV433.....	9	CEV540-5.....	23	CEV634-1A.....	49
CEV433-1.....	8	CEV540-10.....	24	CEV635-1.....	32
CEV433DU.....	6	CEV544.....	23	CEV636-1A.....	49
CEV433H.....	8	CEV544IU.....	22	CEV637I.....	49
CEV435-1.....	33	CEV560.....	38	CEV638T5NU.....	46
CEV435-2.....	33	CEV560-1.....	38	CEV640-10.....	24
CEV438.....	59	CEV565.....	38	CEV641DU.....	17
CEV438B.....	59	CEV566.....	38	CEV641FU.....	17
CEV438D5.....	59	CEV567.....	38	CEV641T1U.....	18
CEV440.....	60	CEV568.....	38	CEV642T1U.....	18
CEV440P6.....	60	CEV570.....	31	CEV642TFU.....	17
CEV442.....	25	CEV580-2.....	50	CEV643IU.....	17
CEV442E.....	25	CEV580-3.....	52	CEV646B5U.....	43, 48
CEV442Z.....	25	CEV601-1.....	20	CEV647B5U.....	43, 48
CEV443.....	54	CEV601P.....	20	CEV648B5U.....	43, 48
CEV444.....	69	CEV602-1.....	10	CEV649-10B.....	45
CEV445-5.....	25	CEV603I.....	39	CEV649B5U.....	43, 48
CEV445IU.....	25	CEV604T1U.....	26	CEV649BGPU.....	27
CEV448.....	60	CEV604T3U.....	26	CEV649C5U.....	27, 43, 48
CEV449.....	60	CEV605BGPU.....	27	CEV649C5.....	7
CEV449G.....	60	CEV605GPU.....	27	CEV649GPU.....	27
CEV453.....	55	CEV605GTU.....	26	CEV649M5.....	7
CEV453A.....	55	CEV605GU.....	26	CEV649P5.....	47
CEV455.....	60	CEV605T1U.....	26	CEV649T5BU.....	43, 48
CEV456.....	59	CEV605T3U.....	26	CEV649T5NU.....	46
CEV460U.....	5	CEV605T5U.....	26	CEV649X5.....	45
CEV470T1U.....	5	CEV606-1.....	14	CEV650T5D.....	23
CEV470U.....	3	CEV607CU.....	27	CEV650TU.....	22
CEV480.....	64	CEV607DU.....	27	CEV651-1.....	20
CEV500-A10.....	24	CEV607T1U.....	27	CEV652I.....	20
CEV500T5.....	23	CEV607U.....	27	CEV653-1.....	29
CEV500TU.....	22	CEV608-1.....	34	CEV658-10.....	44
CEV501.....	21	CEV608D.....	34	CEV660T1U.....	4
CEV501-2.....	21	CEV610-1.....	8	CEV661IU.....	4
CEV502.....	11	CEV611T1U.....	14	CEV662IU.....	4
CEV503.....	29	CEV612DU.....	16	CEV669B.....	49, 51
CEV504PS.....	85	CEV612T1U.....	14	CEV669BG.....	49, 51
CEV504S3.....	85	CEV613T1U.....	14	CEV670DU.....	6
CEV504S5.....	85	CEV614T1U.....	12	CEV670T1U.....	2
CEV505.....	33	CEV615T1U.....	12	CEV671T1U.....	4
CEV506.....	15	CEV616IU.....	16	CEV700T5.....	23
CEV507.....	9	CEV617IU.....	12	CEV700T5D.....	23
CEV508.....	35	CEV620T1U.....	12	CEV720.....	47
CEV511.....	15	CEV621T1U.....	10	CEV720B.....	47
CEV512.....	15	CEV622IU.....	39	CEV720BRUC.....	42
CEV514MU.....	28	CEV623T1U.....	26	CEV720RU.....	42
CEV515.....	13	CEV624IU.....	12	CEV721.....	11
CEV521.....	11	CEV625CU.....	2	CEV725.....	3

CEV728BR.....	44	CEV6795B.....	49, 51	CEV8649T5NU.....	46
CEV728BR5U.....	42	CEV7145A3.....	63	CEV8649T5TBU.....	43, 48
CEV728T5.....	41	CEV7145A5.....	63	CEV8649X5.....	45
CEV728T10.....	44	CEV7145A6.....	63	CEV8650TU.....	22
CEV738-5.....	23	CEV7145A7.....	63	CEV8660T1U.....	4
CEV738D5.....	23	CEV7145B6.....	63	CEV8670DU.....	6
CEV738T5.....	23	CEV7145B7.....	63	CEV8826L.....	62, 76
CEV740-5.....	23	CEV7285RU.....	42	CEV8826M.....	62, 76
CEV744T5.....	23	CEV7625I.....	7	CEV8826P.....	62, 76
CEV760.....	5	CEV7660I.....	7	CEV8826V.....	62, 76
CEV804.....	28	CEV8103D5.....	58	CEV9103D5.....	58
CEV805I.....	66	CEV8103G3.....	77	CEV9103G5.....	58
CEV815.....	3	CEV8103G5.....	58	CEV9104IU.....	26
CEV821G.....	11	CEV8104IU.....	26	CEV9107IU.....	69
CEV825.....	53	CEV8107IU.....	69	CEV9134BG.....	50
CEV830.....	52	CEV8126P.....	63	CEV9145A3.....	63
CEV834.....	50	CEV8126PM.....	63	CEV9145A5.....	63
CEV890.....	77	CEV8126V.....	63	CEV9145A6.....	63
CEV890C.....	77	CEV8143A3.....	61, 75	CEV9145A7.....	63
CEV891IU.....	79	CEV8143A4.....	61, 75	CEV9145B6.....	63
CEV891MU.....	80	CEV8143A7.....	61, 75	CEV9145B7.....	63
CEV892.....	81	CEV8143B7.....	61, 75	CEV9150A3.....	64
CEV893IU.....	79	CEV8145A3.....	63	CEV9150A4.....	64
CEV894IU.....	79	CEV8145A5.....	63	CEV9150A6.....	64
CEV894MU.....	80	CEV8145A6.....	63	CEV9150A7.....	64
CEV895IU.....	79	CEV8145A7.....	63	CEV9150B6.....	64
CEV895MU.....	80	CEV8145B6.....	63	CEV9150B7.....	64
CEV896.....	77	CEV8145B7.....	63	CEV9172A3.....	65
CEV896IU.....	79	CEV8407IU.....	4	CEV9172A4.....	65
CEV897.....	78	CEV8434I.....	8	CEV9172A6.....	65
CEV897BG.....	78	CEV8445IU.....	25	CEV9172A7.....	65
CEV897C.....	78	CEV8538TU.....	22	CEV9172B6.....	65
CEV898.....	77	CEV8544IU.....	22	CEV9172B7.....	65
CEV898IU.....	79	CEV8604T1U.....	26	CEV9229-1.....	53
CEV898L.....	77	CEV8605GTU.....	26	CEV9407IU.....	4
CEV899IU.....	79	CEV8605GU.....	26	CEV9445IU.....	25
CEV899U.....	80	CEV8605T1U.....	26	CEV9500TU.....	22
CEV904.....	28	CEV8607T1U.....	27	CEV9544IU.....	22
CEV915G.....	3	CEV8611T1U.....	14	CEV9604T1U.....	26
CEV965.....	23	CEV8612DU.....	16	CEV9605GTU.....	26
CEV980-3.....	56	CEV8613T1U.....	14	CEV9605GU.....	26
CEV980A.....	56	CEV8615T1U.....	12	CEV9605T1U.....	26
CEV980B.....	56	CEV8616IU.....	16	CEV9607T1U.....	27
CEV984.....	56	CEV8621T1U.....	10	CEV9611T1U.....	14
CEV985.....	56	CEV8623T1U.....	26	CEV9612DU.....	16
CEV992.....	81	CEV8624IU.....	12	CEV9614T1U.....	12
CEV992B.....	72	CEV8625CU.....	2	CEV9615T1U.....	12
CEV992C.....	72	CEV8625DU.....	2	CEV9617IU.....	12
CEV992D.....	81	CEV8625EU.....	6	CEV9620T1U.....	12
CEV995MU.....	81	CEV8625T1U.....	2	CEV9623T1U.....	26
CEV996MU.....	81	CEV8626I.....	30	CEV9624IU.....	12
CEV997MU.....	81	CEV8627I.....	30	CEV9625CU.....	2
CEV998MU.....	81	CEV8641T1U.....	18	CEV9625DU.....	2
CEV1019T5BU.....	40	CEV8642T1U.....	18	CEV9625T1U.....	2
CEV1019T5NUC.....	46	CEV8646B5U.....	43, 48	CEV9627I.....	30
CEV1039G5U.....	40	CEV8647B5U.....	43, 48	CEV9630-1A.....	51
CEV1039T5NU.....	46	CEV8648B5U.....	43, 48	CEV9634-1A.....	49
CEV1039T5TBU.....	40	CEV8649B5U.....	43, 48	CEV9641T1U.....	18
CEV1059-10B.....	44	CEV8649BGPU.....	27	CEV9642T1U.....	18
CEV2295A.....	53	CEV8649C5U.....	27, 43, 48	CEV9646B5U.....	43, 48
CEV6285R.....	41	CEV8649GPU.....	27	CEV9647B5U.....	43, 48
CEV6385RU.....	41	CEV8649P5.....	47	CEV9648B5U.....	43, 48





Integra MicroFrance® product reference number	Manufacturer part number equivalent
<b>CM217S</b>	<b>CM217S</b>
<b>MCO731</b>	<b>430490</b>
<b>SM63</b>	<b>SM63</b>
<b>SM64</b>	<b>SM64</b>
<b>SM65</b>	<b>SM65</b>



## Manufacturers :



Integra MicroFrance SAS (France)  
Le Pavillon ■ 03160 Saint Aubin le Monial ■ France



HEMODIA SAS  
85, rue du Chêne Vert  
31670 LABEGE ■ France



Swann Morton Ltd  
Penn Works ■ Owlerton Green ■ Sheffield  
South Yorkshire S6 2BJ ■ Royaume Uni

## Manufacturers with CE class I products only :



Chemische Fabrik Dr. Weigert GmbH & Co.  
Mühlenhagen  
D-20539 Hamburg ■ Allemagne

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Consult product labels and inserts for any indications, contraindications, hazards, warnings, precautions, and instructions for use.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

### Additional information for EMEA Customers only :

Products mentioned in this document are CE class I, IIa or IIb devices. Contact Integra should you need any additional information on devices classification. All the medical devices mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives.

References that are not medical devices are identified with the indication "NOT A MEDICAL DEVICE".

Products marked with the symbol <sup>®</sup> are distributed by Integra MicroFrance SAS (France)

Date of last révision : 12/2020

## For more information, please contact customer service :

### Europe, Middle-East, Africa

International: +33 (0) 437 47 59 50 ■ +33 (0) 437 47 59 25 (Fax) ■ [csemea@integralife.com](mailto:csemea@integralife.com)  
France: +33 (0) 437 47 59 10 ■ +33 (0) 437 47 59 29 (Fax) ■ [custservfrance@integralife.com](mailto:custservfrance@integralife.com)  
United Kingdom: +44 (0) 1264 312 725 ■ +44 (0) 1264 312 821 (Fax) ■ [custsvcs.uk@integralife.com](mailto:custsvcs.uk@integralife.com)  
Belgium & Luxembourg: +32 (0)2 257 4130 ■ +32 (0)2 253 2466 (Fax) ■ [custsvcbenelux@integralife.com](mailto:custsvcbenelux@integralife.com)  
Netherlands: +31(0) 852 083 167 ■ +31(0) 207 093 627 (Fax) ■ [custsvcnetherlands@integralife.com](mailto:custsvcnetherlands@integralife.com)  
Italy: +39 (0)2 577 8921 ■ +39 (0)2 575 113 71 (Fax) ■ [custsvcitaly@integralife.com](mailto:custsvcitaly@integralife.com)  
[integralife.eu/contact](http://integralife.eu/contact)

Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries.

MicroFrance is a trademark of Integra LifeSciences Corporation or its subsidiaries.

All other trademarks and product names are the property of their respective owners.

© 2020 Integra LifeSciences Corporation. All rights reserved.

0328572-4-EN

DOCUMENT FOR USE IN EUROPE, MIDDLE-EAST AND AFRICA ONLY.